

CnT-ABM10

[Order Now \(/products/ordering/\)](/products/ordering/)

Pen / Strep / Amphotericin B Solution, Single Aliquots

250 x concentrate, sterile solution, ready-to-use single aliquots of penicillin G sodium salt, streptomycin sulfate and amphotericin B (Fungizone). Working concentration is 100 U / mL penicillin, 100 µg / mL streptomycin and 0.25 µg / mL amphotericin B.

TISSUE TYPE	Antibiotics
PACK SIZE	10 x 2.0 mL of frozen concentrate (250 x).
DEFINED	Yes
ACF	Yes
STORAGE / SHELF LIFE	Store at -20°C. See label for expiry date.
SHIPPING CONDITION	CnT-ABM10 is shipped frozen.
COMMENTS	Working concentration is the generally accepted concentration for this antibiotic / antimycotic combination. Some cells may exhibit various degrees of sensitivity to antibiotics or antimycotics.
INTENDED USE	For research use only. Not for use in therapy or diagnostics.
LAST UPDATE	2017-09-27