

## CnT-ST-100

[Order Now \(/products/ordering/\)](/products/ordering/)

### Stain Kit, including Solution I and Solution II

The stain kit can be used to stain fixed cells which were growing on different surfaces. The stain kit consists of two different ready-to-use components. Both stain solutions can be re-used up to five times.

<b>TISSUE TYPE</b>	3D
<b>PACK SIZE</b>	100 mL staining kit
<b>DEFINED</b>	No
<b>ACF</b>	No
<b>STORAGE / SHELF LIFE</b>	Store at room temperature. See labels for expiry date.
<b>QUALITY CONTROL</b>	Stain kit is tested by staining primary human epithelial cells on cell culture plastics and polycarbonate cell culture inserts.
<b>SHIPPING CONDITION</b>	Shipped at ambient temperature.
<b>INTENDED USE</b>	For research use only. Not for use in therapy or diagnostics.
<b>FORM</b>	CnT-ST-100 is a 100 mL staining kit consisting of: <ul style="list-style-type: none"><li>• 1x CnT-ST-I (50 mL, red color, will change to orange during use).</li><li>• 1x CnT-ST-II (50 mL, blue color).</li></ul>
<b>LAST UPDATE</b>	2014-03-21